

## Product datasheet for **RR213114L3V**

### Lsr (NM\_032616) Rat Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Lsr (NM_032616) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Lsr
Synonyms:	Lisch7
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_032616
ORF Size:	1779 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR213114).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_032616.1</a> , <a href="#">NP_116005.1</a>
RefSeq Size:	2097 bp
RefSeq ORF:	1782 bp
Locus ID:	64355
UniProt ID:	<a href="#">Q9WU74</a>
Cytogenetics:	1q21
Gene Summary:	binds apoB, apoE, and triglyceride-rich lipoprotein (chylomicrons and very low density lipoprotein (VLDL)); may mediate the partitioning of dietary lipid between different tissues [RGD, Feb 2006]



[View online »](#)