

Product datasheet for **RR213147L4V**

Cyp2e1 (NM_031543) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Cyp2e1 (NM_031543) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Cyp2e1
Synonyms:	Cyp2e
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_031543
ORF Size:	1479 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR213147).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_031543.1 , NP_113731.1
RefSeq Size:	1624 bp
RefSeq ORF:	1482 bp
Locus ID:	25086
UniProt ID:	P05182
Cytogenetics:	1q41
Gene Summary:	facilitates hepatocyte sensitization to tumor necrosis factor-alpha (TNF-alpha)-mediated cell death; may play a role in xenobiotic metabolism [RGD, Feb 2006]



[View online »](#)