

## Product datasheet for RR213157L4

### Hbegf (NM\_012945) Rat Tagged Lenti ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Hbegf (NM_012945) Rat Tagged Lenti ORF Clone
Tag:	mGFP
Symbol:	Hbegf
Synonyms:	Dtr; GFHB; Hb-egf; Hegfl
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
E. coli Selection:	Chloramphenicol (34 ug/mL)
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR213157).
Restriction Sites:	SgfI-MluI
Cloning Scheme:	

Cloning sites used for ORF Shuttling:



\* The last codon before the Stop codon of the ORF.



[View online »](#)



Locus ID: 25433

UniProt ID: [Q06175](#)

Cytogenetics: 18p11

Gene Summary: activates epidermal growth factor receptor mediated signaling; contributes to neuronal survival; may play a role in neurogenesis during ventral midbrain development [RGD, Feb 2006]