

Product datasheet for **RR213175L4V**

Nhlrc1 (NM_199236) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Nhlrc1 (NM_199236) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Nhlrc1
Synonyms:	Epm2b
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_199236
ORF Size:	1188 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR213175).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_199236.1 , NP_954706.1
RefSeq Size:	1191 bp
RefSeq ORF:	1191 bp
Locus ID:	364682
UniProt ID:	Q6IMG5
Cytogenetics:	17p14
Gene Summary:	human homolog is a putative E3 ubiquitin ligase associated with Lafora progressive myoclonus epilepsy and possibly involved in clearance of polyglucosans from dendrites [RGD, Feb 2006]



[View online »](#)