

Product datasheet for **RR213180**

C1qb (NM_019262) Rat Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: C1qb (NM_019262) Rat Tagged ORF Clone
Tag: Myc-DDK
Symbol: C1qb
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >RR213180 representing NM_019262
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGAAGACACAGTGGAGTGAGATCTTGACACCCCTGTTGCTGCTGCTCCTGGGTTTGCTCCATGTCTCCT
GGGCCAAAGCAGCTGTACTGGGTCCCCTGGCATCCCTGGGGTCCCTGGCATCCCTGGGGTCCCTGGCTC
TGATGGCAAACAGGCACTCCAGGGATAAAAGGAGAGAAAGGGCTCCCCGACTGGCTGGAGACCATGGT
GAGTTAGGAGAGAAAGGGATGCAGGGATCCCTGGGATCCAGGCAAAGTTGGCCCCAAGGTCCCGTCG
GCCCTAAGGGTGCTCCAGGCCCCCTGGACCCCGGTCCTCCAAAGGTGGCTCTGGAGACTACAAGGTAC
CCAGAAAGTAGCCTTCTCTGCCCTGAGGACGGTCAACAGCGCCCTGCGACCAAACAGGCCATTCGCTTC
GAAAAGGTGATACCAATGTTAATGATAACTACGAGCCGCGCAGTGGCAAGTTCACCTGCAAGGTACCTG
GCCTCTACTACTCACCTACCACGCCAGTTCGCGGGAATCTGTGTGTGAACATCGTGCAGCGCCGCGA
CCGAGACCGCATGCAGAAAGTTCTACCTTCTGCGACTATGCCAAAACACCTTCCAGGTACCACGGGT
GGGGTAGTCTTGAAGCTGGAGCAGGAAGAGGTTGTTACCTGCAGGCCACAGACAAGAATCCCTGCTGG
GCGTCGAGGGAGCAATAGCATCTTCACTGGCTTCTGCTTTCCCTGACATGGATGTA

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RR213180 representing NM_019262
Red=Cloning site Green=Tags(s)

MKTQWSEILTPLLLLLLGLLHVSWAQSSCTGSPGIPGVPIPGVPGSDGKPGTPGIKGEKGLPGLAGDHG
ELGEGDAGIPGIPGKVGPKGPVGPKGAPPPGPRGPKGGSDYKATQKVAFSALRVTNSALRPNQAIRF
EKVITNVNDNYEPRSGKFTCKVPLYYFTYHASSRGNLCVNIIVRGRDRDRMQKVLTFCDYAQNTFQVTTG
GVVLKLEQEEVHLQATDKNSLLGVEGANSIFTGFLFPDMDV


TRTRPLEQKLISEEDLAANDILDYKDDDDKV




[View online »](#)

Restriction Sites: SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_019262

ORF Size: 759 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_019262.1</u> , <u>NP_062135.1</u>
RefSeq Size:	1136 bp
RefSeq ORF:	762 bp
Locus ID:	29687
UniProt ID:	<u>P31721</u>
Cytogenetics:	5q36
MW:	26.6 kDa
Gene Summary:	complement protein involved in innate immune response [RGD, Feb 2006]