

Product datasheet for **RR213196L4V**

Sort1 (NM_031767) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Sort1 (NM_031767) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Sort1
Synonyms:	Nt3; Nts3
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_031767
ORF Size:	2475 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR213196).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_031767.1 , NP_113955.1
RefSeq Size:	6690 bp
RefSeq ORF:	2478 bp
Locus ID:	83576
UniProt ID:	O54861
Cytogenetics:	2q34
Gene Summary:	may be involved in sorting luminal proteins in the trans-Golgi [RGD, Feb 2006]



[View online »](#)