

Product datasheet for **RR213219L3V**

Cdc20 (NM_171993) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Cdc20 (NM_171993) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Cdc20
Synonyms:	p55cdc
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_171993
ORF Size:	1497 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR213219).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_171993.1 , NP_741990.1
RefSeq Size:	1777 bp
RefSeq ORF:	1500 bp
Locus ID:	64515
UniProt ID:	Q62623
Cytogenetics:	5q36
Gene Summary:	may play a role in cell cycle regulated protein degradation [RGD, Feb 2006]



[View online »](#)