

## Product datasheet for RR213225L4

### Hagh (NM\_033349) Rat Tagged Lenti ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Hagh (NM_033349) Rat Tagged Lenti ORF Clone
Tag:	mGFP
Symbol:	Hagh
Synonyms:	Glo2; RSP29
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
E. coli Selection:	Chloramphenicol (34 ug/mL)
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR213225).
Restriction Sites:	SgfI-MluI
Cloning Scheme:	

Cloning sites used for ORF Shuttling:



\* The last codon before the Stop codon of the ORF.



[View online »](#)



Locus ID: 24439  
UniProt ID: [O35952](#)  
Cytogenetics: 10q12  
Gene Summary: thiolesterase that catalyzes the hydrolysis of S-D-lactoyl-glutathione; may be involved in the regulation of spermatogenesis [RGD, Feb 2006]