

## Product datasheet for **RR213269L3V**

### Ube2d4 (NM\_031001) Rat Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Ube2d4 (NM_031001) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Ube2d4
Synonyms:	RGD69425; Ube2d2; Ube2d2b
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_031001
ORF Size:	441 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR213269).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_031001.2</a> , <a href="#">NP_112263.1</a>
RefSeq Size:	1499 bp
RefSeq ORF:	444 bp
Locus ID:	79435
UniProt ID:	<a href="#">P70711</a>
Cytogenetics:	14p22
Gene Summary:	ubiquitin-conjugating enzyme (E2 class); mediates the targeting of abnormal and short-lived proteins for degradation [RGD, Feb 2006]



[View online »](#)