

## Product datasheet for **RR213276**

### Cxcl10 (NM\_139089) Rat Tagged ORF Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Cxcl10 (NM\_139089) Rat Tagged ORF Clone  
**Tag:** Myc-DDK  
**Symbol:** Cxcl10  
**Synonyms:** IP-10; Scyb10  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**ORF Nucleotide Sequence:** >RR213276 representing NM\_139089  
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**GCGATCGCC**

ATGAACCCAAGTGTCTGTCGTTCTCTGCCTCGTGTGCTGAGTCTGAGTGGGACTCAAGGGATCCCTC  
TCGCAAGAACGGTGCCTGCACCTGCATCGACTCCATGAACAGACGCTGAGACCCAGGGCCATAGGAAA  
ACTTGAAATCATTCTGCAAGTCTATCCTGTCCGCATGTTGAGATCATTGCCACAATGAAGAAGAACAAT  
GAGAAGAGGTGTCTGAATCCGGAATCTGAGGCCATCAAGAGCTTATTGAAAGCGGTGAGCCAAGAAGGT  
CAAAAAGAGCTCCG

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Protein Sequence:** >RR213276 representing NM\_139089  
Red=Cloning site Green=Tags(s)

MNPSAAVVLCLVLLSLSGTQGIPLARTVRCTCIDFHEQTLRPRAIGKLEIIPASLSCPHEI IATMKKNN  
EKRCLNPESEAIKSLKAVSQRRSKRAP

**TR**TRPLEQKLI SEEDLAANDILDYKDDDDKV

**Restriction Sites:** Sgfl-MluI



[View online »](#)



<b>OTI Annotation:</b>	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u><a href="#">NM_139089.1</a></u> , <u><a href="#">NP_620789.1</a></u>
<b>RefSeq Size:</b>	1133 bp
<b>RefSeq ORF:</b>	297 bp
<b>Locus ID:</b>	245920
<b>UniProt ID:</b>	<u><a href="#">P48973</a></u>
<b>Cytogenetics:</b>	14p22
<b>MW:</b>	10.7 kDa
<b>Gene Summary:</b>	induces DNA synthesis, cell proliferation, and cell migration in vascular smooth muscle cells; may play a role in vascular remodeling [RGD, Feb 2006]