

Product datasheet for RR213277

Gmn (NM_001106112) Rat Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Gmn (NM_001106112) Rat Tagged ORF Clone
Tag: Myc-DDK
Symbol: Gmn
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >RR213277 representing NM_001106112
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCGCGATCGCC

ATGAATCTCAGTATGAAGCAGAAACAGGAAGGAGCCCAAGAGAACGTGAAGAGTAGTCTGTCCCAAGGA
GAACGCTGAAGATGATCCAGCCTTCTGCAGCTGGGTCTCTTGTGGCAGAGAAAAAGATTGCCAAAAGG
CTGTGCCAAAAGGAAACTTTGGGATGACCAATTAACATCTAAGATTTCAAGCTGTGGTCTGAAGCCAAT
GAAAATAACAGTGTGGAGACCTCACCCAGGAAGCATGTGATCTTATGACTAAAGAAAACCCATCTTCTC
AGTATTGGAAAGAAGTGGCAGAACACCGAGGAAAGCTCTCTATGAAGCACTTACGGAGAACGAGAACT
TCATAAAGAAATTGAGCAGAAGGACAGTGAGATTGCCCGCTGAGAAAGGAGAATAAAGACCTGGCAGAA
GTGGCTGAGCACGTACAGTACATGGCAGAGGTAATCGAGAGGCTGAGTAATGAACCTCTGGATAACTTTG
AATCACCGGATAGTCAAGAATTTGATTCTGAAAAAGAACTGTTGAGTATTCTGAACCGGAAGACTCAGA
CACTGGGACATGTCTGAAGAGACTGTCTTCTCTTCCGATGGCACGGCGTGTACA

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RR213277 representing NM_001106112
Red=Cloning site Green=Tags(s)

MNL SMKQKQEGAQENVKSSPVPRRTLKMIQPSAAGSLVGRENELPKGLSKRKLWDDQLTSKISSCGPEAN
ENNSVGDLTQEACDLMTKENPSSQYWKVAEQRRKALYEALTENEKHLKEIEQKDSIARLRKENKDLAE
VAEHVQYMAEVIERLSNEPLDNFESPDSQEFDSEKETVEYSEPEDSDTGTCAEETVSSSSDGTACT

TRTRPLEQKLI SEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI




[View online »](#)

Cloning Scheme:



Plasmid Map:



ACCN: NM_001106112

ORF Size: 618 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001106112.1, NP_001099582.1</u>
RefSeq Size:	1178 bp
RefSeq ORF:	621 bp
Locus ID:	291137
Cytogenetics:	17p11
MW:	23.1 kDa