

## Product datasheet for **RR213278**

### Dapp1 (NM\_001108568) Rat Tagged ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Dapp1 (NM_001108568) Rat Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Dapp1
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin
ORF Nucleotide Sequence:	>RR213278 representing NM_001108568 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGC**C

ATGGGCAGAGCAGAACTTCTAGGAGGGAACATGAGTACCCAGGATCCCTCGGAGCTGTGGGGCCGGGCAG  
ATGGAGGCTCAGACCTGTTCCAGGATTTGGGGTGGTACCATAGCAACCTCACCCGCCATGCTGCAGAGGC  
CCTTCTGCTCTCCAACGGACGCGATGGCAGCTACCTGCTGAGGGACAGCAATGAGCAGACCGGGCTCTAC  
TCTCTCTCCGTGAGAGCCAAAGACTCCGTCAAACACTTTCATGTTGAATACACTGGGTACTCCTTTAAGT  
TTGGCTTTAATGAATATTCATCTCTAAAGGATTTCTGCAAGCATTGCAAAATCAACCTTTGATTGGAAG  
CGAGACAGGGACTCTGATTGTTCTGAAGCATCCCTATCCAAGGCAAGTGGAGGAACCTTCCATTTATGAG  
TCAGTCCGGTCCACACGGCGATGCAGACGGGAAGAACAGAGAATGACCTCGTCCCACTGCACCTTCTC  
TGGGCACCAAAGAAGGCTATCTACCAAGCAGGGAGGTCTGGTAAAGACTTGAAAACGAGATGTTTAC  
TTTGCAAAGGAACGAAGTAAATTTTCAAAGACCAGATGTCTCCAGAACCAATTCGGATCCTAGACCTA  
ACGGAGTGCTCAGCTGTCCAATTTGATTACTCGCAGGATCGCGTAAACTGTTTCTGCTTAGTGTCCCAT  
TCAGGACATTTTATCTCTGTGCAAAGACGGGAGTCAAGCCGATGAATGGATCAAATACTGCGATGGAA  
GTTGTCAAAAATAAGAAAACAGCTGGATCAAGACGACACCATCCGATCTCGGTCAATTTATCTTCAA

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)



<b>OTI Disclaimer:</b>	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
<b>OTI Annotation:</b>	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_001108568.1</a> , <a href="#">NP_001102038.1</a>
<b>RefSeq Size:</b>	3062 bp
<b>RefSeq ORF:</b>	840 bp
<b>Locus ID:</b>	362046
<b>Cytogenetics:</b>	2q44
<b>MW:</b>	32.1 kDa