

Product datasheet for RR213299

Gpr85 (NM_022254) Rat Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Gpr85 (NM_022254) Rat Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Gpr85
Synonyms:	Srep2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>RR213299 representing NM_022254 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGCGAACTATAGCCATGCAGCTGACAACATTTGCAAAATCTCTCGCCTCTAACAGCCTTTCTGAAAC
TGACTTCCTTGGGTTTCATAATAGGAGTCAGTGTGGTGGGCAACCTTCTGATCTCCATTTGCTAGTGAA
AGATAAGACCTTGCATAGAGCTCCTTACTACTCTCTGCTGGATCTGTGCTGCTCAGACATCCTCAGATCT
GCAATTTGTTTTCCATTTGTATTCAACTCTGTCAAAAATGGCTCTACCTGGACTTACGGGACTCTGACTT
GCAAAGTGATTGCCTTTCTGGGGTTTTGTCTGTTTCCACACTGCCTTCATGCTCTTCTGCATCAGCGT
CACCAGATACTTAGCCATCGCCCATCACCGCTTCTATACAAAGAGGCTGACCTTTTGGACGTGTTTGGCT
GTGATCTGCATGGTGTGGACTCTGTCTGTGGCCATGGCATTTCCTCCAGTTTTAGATGTAGGCACCTACT
CATTTCATTAGGGAGGAGGATCAGTGTACCTTCCAACACCGCTCCTTCAGGGCTAACGATTCCCTAGGATT
TATGCTGCTCCTTGTCTCATCCTCCTAGCCACACAGCTTGTCTACCTCAAGCTGATATTTTTGTCCAC
GATCGAAGGAAAATGAAGCCAGTCCAGTTTGTAGCAGCAGTGAGTCAGAACTGGACCTTTCATGGCCCTG
GAGCTAGTGGCCAGGCAGCTGCCAATTGGCTAGCAGGATTTGGAAGGGTCCCACACCACCCACCTTGCT
GGGCATCAGGCAAAATGCGAATACCACAGGCAGAAGACGGCTCTGGTTTTGGATGAGTTCAAAATGGAG
AAAAGAATCAGCAGAATGTTCTATATAATGACTTTCCTTCTAACCTTGTGGGTCCTACCTGGTGG
CCTGCTATTGGAGAGTTTTTGAAGAGGCTGTAGTACCAGGGGATTCTAACAGCCGCTGTCTGGAT
GAGTTTTGCCCAAGCAGGAATCAATCCCTTTGTCTGCATTTTCTCCAACAGGGAGCTGAGGCGCTGTTTC
AGCACAACTTCTTTACTGCAGAAAATCCAGGTTACCAAGGAACTTACTGTGTATA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTAA



[View online »](#)

Protein Sequence: >RR213299 representing NM_022254
Red=Cloning site Green=Tags(s)

MANYSHAADNILQNL SPLTAFLKLTSLGFIIGVSVVGNLLISILLVKDKTLHRAPYYFLDLCCSDILRS
 AICFPFVFN SVKNGSTWYGTLTCKVIAFLGVLSCFHAFMLFCISVTRYLAI AHHRFYTKRLTFWTCLA
 VICMVWTL SVAMAFPPVLDVGTYSFIREEDQCTFQHRSFRANDSLGFMLLLALILLATQLVYLKLIFFVH
 DRRKMKPVQFVAAVSQNWTFHGPASGQAAANWLAGFGRGPTPTLLGIRQNAVNTGRRRLLVLDEFKME
 KRISRMYIMTFLFLTLWGPYLVACYWRVVFARGPVVPGGFLTAAVWMSFAQAGINPFVCIFSNRELRRCF
 STLLLYCRKSRLPREPYCVI

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/ja2092_e02.zip

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_022254

ORF Size: 1110 bp

OTI Disclaimer: Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_022254.3](#), [NP_071590.1](#)

RefSeq Size: 2747 bp

RefSeq ORF: 1113 bp

Locus ID: 64020


UniProt ID: [P60895](#)

Cytogenetics: 4q21

MW: 42.4 kDa

Gene Summary: member of a family of G-protein coupled receptors displaying high conservation across species; rat and human orthologs display complete amino acid sequence identity [RGD, Feb 2006]

Product images:



Circular map for RR213299