

Product datasheet for **RR213341L3V**

Tnni1 (NM_017184) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Tnni1 (NM_017184) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Tnni1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_017184
ORF Size:	561 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR213341).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_017184.1 , NP_058880.1
RefSeq Size:	821 bp
RefSeq ORF:	564 bp
Locus ID:	29388
UniProt ID:	P13413
Cytogenetics:	13q13
Gene Summary:	thin filament contractile regulatory protein; inhibits contractile interaction of myosin and actin [RGD, Feb 2006]



[View online »](#)