

Product datasheet for **RR213385**

Pex3 (NM_031350) Rat Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Pex3 (NM_031350) Rat Tagged ORF Clone
Tag: Myc-DDK
Symbol: Pex3
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >RR213385 representing NM_031350
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGCTGAGATCAATGTGGAATTTCTGAAACGTCACAAAAAGAAATGCATCTTCCTGGGCACAGTCCTTG
GAGGAGTCTATATCCTGGGAAAAATGGACAGAAAAAACTTAGAGAAATACAAGAAAGAGAAGCTGCAGA
GTACATTGCCAAGCTCGGCACAGTACCATTTGAAAGTAACCCAGAGGACTTGCAATATGACAGTGTG
TCCATGCTGCCAACACTGAGAGAGGCCTTAATGCAGCAGCTCAACTCGGAGAGTCTGACAGCTTTGTTGA
AAAACAGGCCTTCAAACAAGCTGGAAATATGGGAGGATCTAAAGATAATAAGTTTACAAGAAGTATTGT
AGCTGTGTACAGCACATGTATGCTGGTGGTTCTTTTGGGAGTCCAGCTGAACATAATTGGTGGATATATT
TACCTGGATAATGCAACAGTTGGCAAAAAATGGCACTAGCATCCTAGCTCCCCAGATGTACAGCAACAGT
ATTTATCAAGTATTCAACACCTCCTTGGAGATGGCCTGACAGAGTTGGTCACTGTCATTAACAAGCTGT
GCAGAGGATCTTAGGAAGTATTTCTCTTAAACATTCCTTGTCCCTTTGGACTTGGAGCAGAAAAGTAAA
GAAATAAGAACTCTTGTGAGCAACATCGATCTTGCTGGAATGATAAAGACGCATCCAAGTCTTCCCTGT
GCCACTACATGATGCCAGACGAGGAGACCGCTGGCAGCGCAGGCTTATGGGCTTCTCCAGAGATAT
TACCACTATTAACCTTCTCAATGAGACAAGAGACATGCTGGAAAGTCCAGATTTCAAGTACAGTTCTGAAT
ACCTGTTTAAACAGAGGCTTCAAGTCCGCTTCTAGACAATATGGCGGAGTTCTTCCGACCCACTGAGCAGG
ACCTGCAGCATGGTAACTCCATAAACAGTCTGTCCAGTGTGAGCCTCCCTTTAGCCAAGATAATCCCAT
AGTGAATGGACAGATCCACTCCGTCTGCAGTGGACACCCAGTCACTTCTGTCAGGATCTGCTGATGATG
GAGCAAGTGAAGACTTTGCTGCCAACGTGTACGAAGCATTACAGACCCCCCAACAAGTGGAGAAA

ACGCGTACGCGGCGGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)


Protein Sequence: >RR213385 representing NM_031350
Red=Cloning site Green=Tags(s)

MLRSMWNFLKRHKKCI FLGTVLGGVYILGKYGQKLR EIQERAAEYIAQARRQYHFESNQR TCNMTVL
 SMLPTLREALMQLNSESL TALLKNRPSNKLEIWEDLKIISFTRSI VAVYSTCMLVLLRVQLNIIGGYI
 YLDNATVGKNGT SILAPDPVQQYLSSIQHLLGDGL TELVTVIKQAVQRILGSISLKHSLSLLDLEQK LK
 EIRTLVEQHRSCWNDKASKSSLCHYMPDEETPLAAQAYGLSPRDITTIKLLNETRDMLESPDFSTVLN
 TCLNRGFSRLLDNMAEFFRPTEQDLQHNSINLSSVSLPLAKIIPVINGQIHSVCSETPSHFVQDLLMM
 EQVKDFAANVYEAFSTPQQLEK


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_031350

ORF Size:	1116 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_031350.1 , NP_112640.1
RefSeq Size:	1952 bp
RefSeq ORF:	1119 bp
Locus ID:	83519
UniProt ID:	Q9JJK4
Cytogenetics:	1p13
MW:	42.2 kDa
Gene Summary:	integral membrane protein that binds Pex19p; plays a role in assembly of peroxisomal membrane vesicles [RGD, Feb 2006]