

Product datasheet for **RR213411L3V**

Psmb10 (NM_001025637) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Psmb10 (NM_001025637) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Psmb10
Synonyms:	MGC112890
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001025637
ORF Size:	819 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR213411).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001025637.1 , NP_001020808.1
RefSeq Size:	1021 bp
RefSeq ORF:	822 bp
Locus ID:	291983
UniProt ID:	Q4KM35
Cytogenetics:	19q12



[View online »](#)

Gene Summary:

The proteasome is a multicatalytic proteinase complex which is characterized by its ability to cleave peptides with Arg, Phe, Tyr, Leu, and Glu adjacent to the leaving group at neutral or slightly basic pH. The proteasome has an ATP-dependent proteolytic activity. This subunit is involved in antigen processing to generate class I binding peptides (By similarity).
[UniProtKB/Swiss-Prot Function]