

Product datasheet for **RR213416L4V**

Nr1i2 (NM_052980) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Nr1i2 (NM_052980) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Nr1i2
Synonyms:	PXR
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_052980
ORF Size:	1293 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR213416).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_052980.2
RefSeq Size:	2360 bp
RefSeq ORF:	1296 bp
Locus ID:	84385
UniProt ID:	Q9R1A7
Cytogenetics:	11q21
Gene Summary:	involved in the induction of cytochrome P4503A (CYP3A) [RGD, Feb 2006]



[View online »](#)