

Product datasheet for **RR213437L3V**

Gdf6 (NM_001013038) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Gdf6 (NM_001013038) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Gdf6
Synonyms:	BMP-13; gdf16
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001013038
ORF Size:	1356 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR213437).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001013038.1 , NP_001013056.1
RefSeq Size:	1576 bp
RefSeq ORF:	1359 bp
Locus ID:	252834
UniProt ID:	Q6HA10
Cytogenetics:	5q13
Gene Summary:	member of the bone morphogenetic protein (BMP) family and the TGF-beta superfamily; regulates cell growth and differentiation; involved with tendon and ligament formation [RGD, Feb 2006]



[View online »](#)