

Product datasheet for **RR213443L3V**

Nr4a2 (NM_019328) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Nr4a2 (NM_019328) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Nr4a2
Synonyms:	Nurr1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_019328
ORF Size:	1794 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR213443).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_019328.3 , NP_062201.2
RefSeq Size:	2582 bp
RefSeq ORF:	1797 bp
Locus ID:	54278
UniProt ID:	Q07917
Cytogenetics:	3q21
Gene Summary:	nuclear receptor and transcription factor; plays a role in development and maintenance of neurons synthesizing the neurotransmitter dopamine [RGD, Feb 2006]



[View online »](#)