

Product datasheet for **RR213468**

Slc5a6 (NM_130746) Rat Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Slc5a6 (NM_130746) Rat Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Slc5a6
Synonyms:	SMVT
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin



[View online »](#)

ORF Nucleotide
Sequence:

>RR213468 representing NM_130746
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCGCGATCGCC

ATGACTGTGGCGAGCACTGCGGCTCCCTCCTATAACAACCTCAGACACCAACCGGGTTATTTCCACCTTTT
CTGTTGTGGACTATGTGGTGTGGCCTGTTGCTGGTTCTCTCCCTTGTCATTGGGCTCTATCATGCTTG
CCGAGGCTGGGGCCGCCATACTGTTGGTGGAGCTGCTGATGGCAGACCGAAAAATGGGCTGCCTTCTGTG
GCACTGTCCTTGTGGCCACCTTCCAGTCAGCCGTGGCGATCCTTGGGGGACCAGCTGAGATCTACCGAT
TTGGAACCCAGTATTGGTTCCTGGGGTCTCCTACTTCTGGGGCTTCTGATCCCTGCTCACATCTTCAT
CCCTGTCTTACCGCCTGCATCTTACCAGCGCCTATGAGTACCTAGAGCTCCGCTTCAATAAAGCAGTC
CGGATCTGTGGGACTGTGACCTTCATCTTTCAGATGGTGGTCTACATGGGAGTGGCTCTCTATGCACCGT
CCCTGGCCCTCAATGCAGTGACTGGATTTGATCTATGGCTGTGAGTCTGGCCCTAGGAATTGTCTGCAA
CATCTACACTGCACTGGGTGGGCTGAAGCTGTGATCTGGACAGATGTGTTCCAAACGTTGATCATGTTT
CTAGGGCAGCTGGTGGTTATCATTGTAGGAGCGGCCAAAAGTGGGCGGCTTGGGGCATGTATGGGCTGTGG
CTTCCCAGCATGGCCTCATCTCTGGGATTGAGCTGGATCCTGACCCCTTGTACGTCATACTTTCTGGAC
TTTGGCCTTTGGCGCGCTTTCATGATGCTCTCCTTATATGGAGTGAATCAGGCTCAGGTGCAGCGCTAC
CTCAGTCCCCTCAGAGAAGGCTGCTGTGCTCTCCTGCTATGCAGTGTTCCTTGGCAGCAAGTGGCC
TCTGCATGAGCTGCCTCATTGGTCTGGTGTGTTGCCTATTATAAGAAATATAGTATGAGCCCCAACA
AGAGCAAGCAGCCCTGACCAGTTAGTCTCTATTTGTGATGGACCTCCTGAAAGACATGCCAGGGCTG
CCCGGGCTCTTCGTTGCTTGCCTCTTTCAGTGGATCTCTCAGCACCATATCCTCCGATTCATTCAGTGG
CCACTGTCACAATGGAAGACCTGATCCAGCCCTGGTCCCTCAGTTGACTGAAACCCGGGCCATCATGCT
CTCTCGAAGCCTTGCCTTTGCCTATGGCTAGTTTGCCTGGGAATGGCCTATGTTTCTCTCATCTGGGA
TCAGTGCTCCAGGCAGCACTCAGCATCTTGGCATGGTTGGAGGGCCACTGCTGGGACTCTTCTGCCTTG
GGATGTTCTTCCGTGTGCCAACCTCTTGGTGCCATTGTAGGCTGTTGACGGGACTCACCATGGCCTT
CTGGATTGGCATTGGGAGCATAGTAAGCAGGATGAGTTCTGCTGCAGCATCCCCTCCCCTTAATGGATCC
AGCTCCTTCTGCCAGCAATCTAACTGTCGCCACTGTGACCACCCTGATGCCTTCTACCCTGTCCAAGC
CCACAGGACTGCAGCAGTTCTATCCCTGTCCTATTTATGGTACAGTGCGCACAACTCCACCACAGTCAT
TGCAGTGGCCTGATTGTCAGTCTGCTCACTGGGGGAATGCGGGCCGGTCCCTGAACCTGGCACCATT
TACCCTGTGTGCCAAAGCTCCTTGCACTCCTGCCCTTATCCTGCCAGAAGCGGCTTGTGCTGGAGAAGCC
ACAACCAGGATATCCCTGTGGTCACCAACCTGTTTCCAGAGAAGATGGGGAACGGAGCGCTGCAGGACAG
CAGAGACAAAGAGAGGATGGCTGAAGATGGCCTTGTCATCAGCCATGCAGCCCCACCTACATTGTCCAA
GAGACGTCCTG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RR213468 representing NM_130746
 Red=Cloning site Green=Tags(s)


MTVASTAAPSYYTSDTNRVISTFSVVDYVVFGLLLVLSLVIGLYHACRGWGRHTVGELLMADRKMGLCPV
 ALSLLATFQSAVAILGGPAEIYRFGTQYWFLGCSYFLGLLIPAHIFIPVFYRLHLTSAYEYLELRFNKAV
 RICGTVTFIFQMVVYMGVALYAPSLALNAVTFGDLWLSVLALGIVCNIYALGGLKAVIWDVVFQTLIMF
 LGQLVVIIVGAAKVGGGLGHVWAVASQHGLISGIELDPDPFVRHTFWTLAFGGVFMMLSLYGVNQAQVQRY
 LSSHSEKAAVLSCYAVFPCQQVALCMSCLIGLVMFAYYKKYSMSPPQEQEAPDQLVLYFVMDLLKDMPLG
 PGLFVAACLFSGSLSTISSAFNSLATVTMEDLIQPWFQLTETRAIMLSRSLAFAYGLVCLGMAYVSSHG
 SVLQAALSIFGMVGGPLLGLFCLGMFFPCANPLGAIVGLLTLGLTMAFWIGIGSIVSRMSSAAASPPLNGS
 SSFLPSNLTVATVTTLPSTLSKPTGLQQFYLSYLWYSAHNSTTVIAGLIVSLLTGGMRGRSLNPGTI
 YPVLPKLLALLPLSCQKRLCWRSHNQDIPVVTNLFPEKMGNGALQDSRDKERMAEDGLVHQPCSPYIVQ
 ETSL

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Restriction Sites:

SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_130746

ORF Size: 1902 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_130746.1](#), [NP_570102.1](#)

RefSeq Size: 3091 bp

RefSeq ORF: 1905 bp

Locus ID: 170551

UniProt ID: [O70247](#)

Cytogenetics: 6q14

MW: 68.6 kDa

Gene Summary: mediates sodium-dependent transport of pantothenate, biotin, and lipoate; plays a role in intestinal uptake of biotin and in transplacental vitamin transport [RGD, Feb 2006]