

Product datasheet for **RR213488**

Pmel (NM_001108729) Rat Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Pmel (NM_001108729) Rat Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Pmel
Synonyms:	Si; Silv
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin
ORF Nucleotide Sequence:	>RR213488 representing NM_001108729 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGGTGTCCAGAGAAGGTGCTTCCTCCTGTGCTTGTGCTGGGTGCTCTGCTGGCTTTGGGGTCTATAG
AAGGATCCAGGAATCAGAACTGGCATGGCGTCTCGAGGCAACTTGAACAAAGTCTGGAATAAGCAGCT
GTACCCCTGAGTGGACAGAGGTGCAGGGTCTAACTGCTGGAGAGGTGGCCAGGTGTCTCTGAAGGTCAGG
AATGATGGCCCTACACTGGTTGGTGAAATACCTCCTTTCCATTGCTCTGCACTTCCCTGGAAGTCAAA
AGGTAAGTCCCGGATGGTCAAGTTATCTGGTCAACAACACCATCATCAATGGGAGCCAGGTGTGGGGAGG
ACAGCCAGTGTATCCACGGGAACCTGATGATGCCTGCATTTTCCCTGACGGTGGACCCTGCCCATCTGGT
CCCAAGCCTCCGAGGAGAAGCTTTGTTTATGTTTGAAGACCTGGGGACAATACTGGCAAGTTCTGGGGG
GCCCGGAGTCCAAGCTGAGCATTCCAACGGGTACGCAAGGCTGGGCACACACACCATGGAAGTGACTGT
CTACCACCGACGGGGTCCCAGAGCTACGTGCCCTGGCTCACTCCAGTTCACCTTCCACCATCACTGGG
TCCGTTAGCCGTCTACTGGATGACACAGACACCAATAATGCTTGTGAAGAGACAAGTCCCCTGGATTGTG
TTCTATATCGATATGGTCTTTTCCCTCACCCCTGGATATCGTCCAGGGATTGAAAGTCTGAGATCCT
GCAGGCTGTGCCATCCAGTGAGGGGGACGCCCTTGGAGCTGACTGTGCTCCTGCCGAGGCGGGCTACCCAAG
GAAGCCTGCATGGATATTTTCATCACAGGGTGCCAGCCCTGCCAGAGGCTGTGCCAGCCTGTGCCAC
CAAGCCCAGACTGCCAGCTGTTCTACACCAAATACTGAAAGGCGGCTTGGGGACTTATTGCCTCAATGT
GTCTTTGGCTGATGCCAACAGCCTGGCAGTGGCCAGCACCCAACCTGTTGTTCTGGTCAAGAAGGCAGC
CTTGGGCAGGCTCCCTTGTAGGTGTCTTACTGGTGTGGTGGCTGTGGTCTTGCATCTCTGATAT
ATAGGCACAGACTTAAGAAGCAGGACTCAGTTTCCCAAACGCCACATGGCAGCACTCACTGGCTGCGCCT
GCCTCCAGTCTTCTGCGCTCGCCGCTCGGTGAAAGCAGTCCGCTCCTCAGCGGGCAGCAGTT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >RR213488 representing NM_001108729
 Red=Cloning site Green=Tags(s)


MGVQRRFCFLPVLVLGALLALGSIEGSRNQNHGVSRLVTKVWNKQLYPEWTEVQGSNCWRGGQVSLKVR
 NDGPTLVGANTSFSIALHFPGSQKVLDPGQVIWVNTIINGSQVWGGQPYPREPDDACIFPDGGPCPSG
 PKPPRRSFVYVWKTWQYVQVLGGPESKLSIPTGHARLGTHTMEVTVYHRRGSQSYVPLAHSSTFTITG
 SVSRLDDTDTIMLVKRQVPLDCVLYRYSFSLTLDIVQGIESAELQAVPSSEGDFAFELTVSCRGLPK
 EACMDISSPGCQPPAQRLLCQVPPSPDCQLVLHQILKGLGTYCLNVSLADANSLAVASTQLVVPQGEGS
 LGQAPLLVGVLLVLVAVVLASLIYRHLKKQDSVSQTPHGSTHWLRLPPVFCARRLGESSPLLSGQQV

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Restriction Sites:

SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_001108729

ORF Size:	1254 bp
OTI Disclaimer:	<p>Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.</p> <p>The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info</p>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001108729.1, NP_001102199.1</u>
RefSeq Size:	1666 bp
RefSeq ORF:	1257 bp
Locus ID:	362818
Cytogenetics:	7q11
MW:	45.3 kDa