

## Product datasheet for RR213490L4

### Aadat (NM\_017193) Rat Tagged Lenti ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Aadat (NM_017193) Rat Tagged Lenti ORF Clone
Tag:	mGFP
Symbol:	Aadat
Synonyms:	Kat2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
E. coli Selection:	Chloramphenicol (34 ug/mL)
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR213490).
Restriction Sites:	SgfI-MluI
Cloning Scheme:	

Cloning sites used for ORF Shuttling:



\* The last codon before the Stop codon of the ORF.



[View online »](#)



Locus ID: 29416  
UniProt ID: [Q64602](#)  
Cytogenetics: 16p12  
Gene Summary: endogenous modulator of glutamatergic neurotransmission with kynurenine aminotransferase (KAT) activity [RGD, Feb 2006]