

Product datasheet for **RR213501L4V**

Steap3 (NM_133314) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Steap3 (NM_133314) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Steap3
Synonyms:	MGC93147
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_133314
ORF Size:	1464 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR213501).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_133314.1 , NP_579848.1
RefSeq Size:	2707 bp
RefSeq ORF:	1467 bp
Locus ID:	170824
UniProt ID:	Q5RKL5
Cytogenetics:	13q11
Gene Summary:	induces apoptosis; inhibits prostate tumor growth in some cell lines [RGD, Feb 2006]



[View online »](#)