

Product datasheet for **RR213535L3V**

Creld1 (NM_001024783) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Creld1 (NM_001024783) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Creld1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001024783
ORF Size:	1260 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR213535).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001024783.1 , NP_001019954.1
RefSeq Size:	1900 bp
RefSeq ORF:	1263 bp
Locus ID:	312638
UniProt ID:	Q4V7F2
Cytogenetics:	4q42
Gene Summary:	Protein disulfide isomerase (By similarity). Promotes the localization of acetylcholine receptors (AChRs) to the plasma membrane (By similarity).[UniProtKB/Swiss-Prot Function]



[View online »](#)