

Product datasheet for **RR213542L4V**

Procr (NM_001025733) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Procr (NM_001025733) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Procr
Synonyms:	MGC114418; MGC116255
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001025733
ORF Size:	723 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR213542).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001025733.2 , NP_001020904.1
RefSeq Size:	1400 bp
RefSeq ORF:	726 bp
Locus ID:	362248
UniProt ID:	Q4V8I1
Cytogenetics:	3q42
Gene Summary:	Binds activated protein C. Enhances protein C activation by the thrombin-thrombomodulin complex; plays a role in the protein C pathway controlling blood coagulation (By similarity). [UniProtKB/Swiss-Prot Function]



[View online »](#)