

Product datasheet for RR213563L4

Fam161a (NM_001013876) Rat Tagged Lenti ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Fam161a (NM_001013876) Rat Tagged Lenti ORF Clone
Tag:	mGFP
Symbol:	Fam161a
Synonyms:	RGD1304999
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
E. coli Selection:	Chloramphenicol (34 ug/mL)
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR213563).
Restriction Sites:	SgfI-MluI
Cloning Scheme:	

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF.

ACCN:	NM_001013876
ORF Size:	1488 bp



[View online »](#)

