

## Product datasheet for **RR213572L3V**

### **Pnlip (NM\_013161) Rat Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	Pnlip (NM_013161) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Pnlip
Synonyms:	PANLI
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_013161
ORF Size:	1395 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR213572).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_013161.2</a> , <a href="#">NP_037293.1</a>
RefSeq Size:	1491 bp
RefSeq ORF:	1398 bp
Locus ID:	25702
UniProt ID:	<a href="#">P27657</a>
Cytogenetics:	1q55



[View online »](#)

**Gene Summary:**

Plays an important role in fat metabolism. It preferentially splits the esters of long-chain fatty acids at positions 1 and 3, producing mainly 2-monoacylglycerol and free fatty acids, and shows considerably higher activity against insoluble emulsified substrates than against soluble ones.[UniProtKB/Swiss-Prot Function]