

## Product datasheet for RR213609L4

### Pik3r2 (NM\_022185) Rat Tagged Lenti ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Pik3r2 (NM_022185) Rat Tagged Lenti ORF Clone
Tag:	mGFP
Symbol:	Pik3r2
Synonyms:	MGC108697
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
E. coli Selection:	Chloramphenicol (34 ug/mL)
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR213609).
Restriction Sites:	SgfI-MluI
Cloning Scheme:	

Cloning sites used for ORF Shuttling:



\* The last codon before the Stop codon of the ORF.

ACCN:	NM_022185
ORF Size:	2166 bp



[View online »](#)

