

Product datasheet for **RR213622L3V**

Pla2g10 (NM_017176) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Pla2g10 (NM_017176) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Pla2g10
Synonyms:	PLA2GX; sPLA2-X
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_017176
ORF Size:	453 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR213622).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_017176.2 , NP_058872.1
RefSeq Size:	910 bp
RefSeq ORF:	456 bp
Locus ID:	29359
UniProt ID:	Q9QZT3
Cytogenetics:	10q11
Gene Summary:	mouse homolog has phospholipase A2 activity [RGD, Feb 2006]



[View online »](#)