

Product datasheet for **RR213634L3V**

St3gal2 (NM_031695) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	St3gal2 (NM_031695) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	St3gal2
Synonyms:	Siat4b; Siat5
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_031695
ORF Size:	1050 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR213634).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_031695.2 , NP_113883.2
RefSeq Size:	3221 bp
RefSeq ORF:	1053 bp
Locus ID:	64442
UniProt ID:	Q11205
Cytogenetics:	19q12
Gene Summary:	Gal beta 1,3GalNAc alpha 2,3-sialyltransferase with potential substrates of asialo-GM1 and GM1 [RGD, Feb 2006]



[View online »](#)