

## Product datasheet for **RR213667L3V**

### Serpine2 (NM\_019197) Rat Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Serpine2 (NM_019197) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Serpine2
Synonyms:	CRG; Gdnpn1; PI-7; Pn-1; Spin4
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_019197
ORF Size:	1191 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR213667).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_019197.1</a> , <a href="#">NP_062070.1</a>
RefSeq Size:	2078 bp
RefSeq ORF:	1194 bp
Locus ID:	29366
UniProt ID:	<a href="#">P07092</a>
Cytogenetics:	9q34
Gene Summary:	displays neurite promoting activity and serine protease inhibitory activity [RGD, Feb 2006]



[View online »](#)