

Product datasheet for **RR213670L3V**

Stk16 (NM_173142) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Stk16 (NM_173142) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Stk16
Synonyms:	F52; PKL12
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_173142
ORF Size:	915 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR213670).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_173142.1 , NP_775165.1
RefSeq Size:	1489 bp
RefSeq ORF:	918 bp
Locus ID:	286927
UniProt ID:	P57760
Cytogenetics:	9q33
Gene Summary:	serine/threonine protein kinase enzyme; may regulate transforming growth factor-beta (Tgfb) and vascular endothelial growth factor transcription [RGD, Feb 2006]



[View online »](#)