

Product datasheet for **RR213678L4V**

II15 (NM_013129) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	II15 (NM_013129) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	II15
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_013129
ORF Size:	486 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR213678).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_013129.2 , NP_037261.1
RefSeq Size:	489 bp
RefSeq ORF:	489 bp
Locus ID:	25670
UniProt ID:	P97604
Cytogenetics:	19q11
Gene Summary:	human homolog is a cytokine involved in the regulation of T and natural killer cell activation and proliferation [RGD, Feb 2006]



[View online »](#)