

Product datasheet for **RR213679**

Hexx1 (NM_001109106) Rat Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Hexx1 (NM_001109106) Rat Tagged ORF Clone
Tag: Myc-DDK
Symbol: Hexx1
Synonyms: RGD1563858
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >RR213679 representing NM_001109106
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGTCTCCAGCCTTCGGGAAGTTGCTCAGCTCCGGGAAAGCAAACCTCACCTGCTCCTTCTCAATCG
AGAGCATTTTAGGACTGGACCAGAAAAAGATTGCGCAACGTAGTAAGACCCACAGACCTGGACAGA
CACCTGCGGCGACTCAGAGAAAGACGGTAACCCACGTCTACATGCCCCAGGTCTCCAGTGAGATTCA
TTTCCTTGTCCAGTGGATCACCCAATGCCTGAAGAAAGGGCTCCCAAATATGAAAATTACTTTTCAGCCT
CAGAAACACACTCTTTGAAAAGAGAGTTGAGTTGGTACAGAGGACGAAGGCCAAGAACCCTTTTACACA
GAACCAGGTGCAAGTACTAGAAAATGTCTTCAAGTGAATGCTATCCTGGCATTGATATCAGAGAGGAC
CTAGCTCAAAAGCTGAATTTAGAGGAGGACAGAAATCCAGATTTGGTTCCAAAACCGTCGAGCAAAGCTGA
AAAGGTCCCGGAGAGAATCACAGTTTCTAATGGCAAAAAGCCCTTCAATCCAGATCTTCTGAAA

ACGCGTACGCGGCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RR213679 representing NM_001109106
Red=Cloning site Green=Tags(s)

MSPSLREVAQLRESKSPCSFSIESILGLDQKKDCATSVRPHRPWTDTCGDSEKDGPNRLHAPGLPSEIS
FPCPVDHPMPEERAPKYENYFSAETHSLKRELSWYRGRPRPTAFTQNQVEVLENVFRMNCYPGIDIRED
LAQKLNLEEDRIQIWFQNRRAKLRSRRESQFLMAKKPFNPDLLK

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Restriction Sites: SgfI-MluI



Cloning Scheme:



Plasmid Map:



ACCN: NM_001109106

ORF Size: 555 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001109106.1</u> , <u>NP_001102576.1</u>
RefSeq Size:	1148 bp
RefSeq ORF:	558 bp
Locus ID:	498575
Cytogenetics:	16p16
MW:	21.6 kDa