

Product datasheet for **RR213696L4V**

Sec22b (NM_001025686) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Sec22b (NM_001025686) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Sec22b
Synonyms:	Sec22l1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001025686
ORF Size:	645 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR213696).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001025686.1 , NP_001020857.1
RefSeq Size:	1465 bp
RefSeq ORF:	648 bp
Locus ID:	310710
UniProt ID:	Q4KM74
Cytogenetics:	2q34
Gene Summary:	SNARE involved in targeting and fusion of ER-derived transport vesicles with the Golgi complex as well as Golgi-derived retrograde transport vesicles with the ER.[UniProtKB/Swiss-Prot Function]



[View online »](#)