

Product datasheet for **RR213699**

Adam30 (NM_001109600) Rat Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Adam30 (NM_001109600) Rat Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Adam30
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin



[View online »](#)


ORF Nucleotide
Sequence:

>RR213699 representing NM_001109600
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGAGGTCAGTGTGGCCTCCCTCCCCAACACCCTGGCTGCTTCTGGCGATGCTTTTCTTGAGGCTA
TGGCCGAAGACCTACTATTTGACCCTGACTGGGGCTTCGACTCGTATGAAATCACCATTCCCAAAGTACT
AAGCTTGAGGAAAGAGGCGCAGGGTGTGACAGTTCTTTGTCATACCTTTACAGATAGAAGCAAGGAG
CGTGTGGTGCCTTGGGCCAAGAAGTTACTGCTGCCAGACACCTGCCAGTTTTCTCCTTACCCAAG
AGGACAGTCTGATAGAGGACTACCCACATATCCAGACCAGTGAATTACGTGGGCTTCGTGGAGGGCTC
TCAGGAAGCCGAATCTACTCTGAGTACATGCATGGGGGTCTGCGAGGCGTCTTGAACATTGACTCCAGA
TACTACAAATCGAACCCCTCAGGGCTCTTCCCGTTTTGAACATGTAGTGTATCTCCTAACTGAGGACA
ACTTCGGCAATGAGACCTGTGGTGTGGTTGATGAAGAAACAGATGGGCAACAGCCGAGGAAGAGGAAGT
GGCTAGAATAAAAACTTTCATCAGTCTACAAGCACAAAGGTAAGTGAATTGGTCATGGTCTTTGAC
CTGAGTAGATATGTGTTACCCGAGGCAATATGACTAGAATGGTAGAGGACGCCATTCTCTTGGCTTCAA
TTGACACACCTACTTTCAAGGATGTCAGCTTGAAGATCAACTGAAGGCGCTTGAAGTCTGGACGAGTT
TGACAAAATACATCTTTATTATCTCTTTATCCCAAGTGTAGGCCAATTTGTGCTTTACAAAAGAGCT
TCGCTGCATCACCGGCTTCTGCGGACTGGGCACATTTATACATTGCAAGGAGATACAGTGTACTG
CGTGGTCTTGGGAAGGGCATGTGAAGAATCCATGCTGGATCCGCAAGCAGCTTCTCAACGTTAATCT
CCTGGGACCTGCCACTGGACTGCTCACGAAATGGGCCATTGTGTGGGATGTATCATGATGAGGAATTC
TGTCAAGTGCAGAGGCAGGAAGAGCTGTGTCATGGGTACAGGACGTACAGGGTTTAGTAACTGTAGTTACC
ATCAGTATTTTTTACACGCGCTTACAAAATGTCTTATGCTACTGACTTCCCAAAGTGCATATAGT
AGAGAGGTGTGGAATAAAATTTGTGAAAACAAGGAGGAGTGTGACTGCGGCTCCGAAAGAAGAGTGCAAA
AGCGACCTGTGCTGCGGGCCAGACTGTAAAGTGAAGAGGGCGTCAACTGCTCCTCCGGGCTCTGTTGTC
ATAAGTCAACTTTCTCCATCGGGCTATGTGTGACAGGAGGAAGAAAATGAATGCGACCTCGCAGAGTA
CTGTGACGGGCTCTCAGGTTTCTGTCCGAAAGACTATTACAAGCAGGATGGGACTCCGTGCAGGTATGGA
GGCTTTTGCTTCCATAAGGGGTGCAGATCCAGTACTTGCAGTGCAGGCCATTTTTGGACGTAGAGCCA
AGGAGGCCCTCGCCAGTGTACGATGCAGTTAACATGATAGGTGACCAGTCTGGGAAGTGTGGCATTGT
AAATGCTACCCGATACACCAAGTGTATAGAGACAACACGATCTGTGGTCCGCTCCAGTGTGTAACATT
AAAGCCCTCTCCAGATTTCTGTATCACACCCTTCTACTGAGCACCCACCTGAAGAAAGACAATCTCAAGT
GCTGGGCACGAGCTACCATGGGACCATGATACCGCACGGGTACCCGACATAGGTGATGTCAACGACGG
TACCTCCTGCGGTCTAACCGGGTGTGCCTAACAGGACTTGCCTTACTACGCCACCCTGAAGTTCGAC
TGTTTGCTAAAAAGTGAATAGCCGGGGTGTGTCACAAACAAAAAGAACTGCCATTGCACATATGGCT
GGGAACCCCATCTGTGTGGAGGTAGGGTACGGTGGGAGCATTGACAGTGGGCCCTGGACCCAGGGC
GAGACTCGAAAGGTTACGTGCTGTACATGCTGCGTACTATTCTGTTTCTTTCTGTGATCGTTGTG
TGTTTTACGCAGCGAATGAAAAATTTGTTTCGTAAGCGCCAGCAAAAAACCCAGCTGAGACTCAGGGGA
AGAAGAGG

AGCGGACCGACGCGTACGCGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCC
TGGATTACAAGGATGACGACGATAAGGTTTAA

Plasmid Map:


ACCN: NM_001109600

ORF Size: 2178 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001109600.1](#), [NP_001103070.1](#)

RefSeq Size: 3441 bp

RefSeq ORF: 2181 bp

Locus ID: 690547
Cytogenetics: 2q34
MW: 82.8 kDa