

Product datasheet for **RR213705**

Cx3cr1 (NM_133534) Rat Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Cx3cr1 (NM_133534) Rat Tagged ORF Clone
Tag: Myc-DDK
Symbol: Cx3cr1
Synonyms: Rbs11
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >RR213705 representing NM_133534
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GCGATCGCC**

ATGCTACCTCCTCCCGAATTGGATCTAGAGAACTTTGAGTATGATGACTCTGCTGAGGCCTGTTATT
TGGGTGACATCGTGGCCTTTGGGACCATCTTCTATCTATTTTCTACTCCCTTGCTTCAGTTCGGTCT
GGTGGGAATCTGTTGGTGGTCTCGCCCTCACCAACAGCCGGAAGTCCAAGAGCATCACTGACATCTAC
CTCCTGAACCTGGCCTTGAGCGACCTGCTCTTTGTGGCCACTTTGCCCTTCTGGACTCACTACCTCATCA
GCCATGAGGGCCTCCACAACGCCATGTGCAAGCTCACGACTGCTTTCTTCTTATTGGCTCTTTGGGGG
CATATTCTTCATCACCGTCATCAGCATCGACCGGTACCTCGCCATCGTCTCGCCGCAACTCCATGAAC
AACCGGACAGTGCAACACGGCGTCACCATCAGTCTGGGCGTCTGGGCGGCGGCCATCTTAGTGGCGTCG
CCCAGTTCATGTTACAAAGAGAAAGGACAACGAATGTTGGGTGATTACCCCGAGGTCCTGCAGGAAAT
CTGGCCCGTGCCTCGCAACTCGGAGGTCAACATCCTGGGCTTCGTCCTGCCCTTGCTTATCATGAGCTTT
TGCTACTCCGCATCGTCCGACGCTGTTTTCTGCAAGAACCGAAGAAGGCCAGAGCCATTAGGCTCA
TCCTCTGGTGGTGTGTTGCTTCTTCTTCTTGACGCCTTACAACATCGTGATTTTCTGGAGACTCT
CAAATTCTACAATTCTTCCCTAGTTGTGGCATGAAGAGGGACCTGAGGTGGGCCCTTAGTGTGACGGAG
ACAGTGGCGTTTAGCCACTGCTGCCTCAACCCCTTATCTACGCTTTCGCTGGGGAAAAGTTCAGAAGGT
ACCTGAGACACCTGTACAACAAGTGCCTGGCCGTCCTGTGCGGTGCTCCTGTCCACGCCGCTTCTCAAC
AGAGTCCCAGAGGAGCAGGCAGGACAGCATTCTGAGCAGCTTGACTCACTACACAAGCGAGGGAGAGGGA
TCTCTCTGCTC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

ORF Size:	1062 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_133534.1 , NP_598218.1
RefSeq Size:	1318 bp
RefSeq ORF:	1065 bp
Locus ID:	171056
UniProt ID:	P35411
Cytogenetics:	8q32
MW:	40.3 kDa
Gene Summary:	member of the rhodopsin family of G-protein coupled receptors; may act as a receptor for a chemokine peptide ligand [RGD, Feb 2006]