

Product datasheet for RR213710

Pdpm (NM_019358) Rat Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Pdpm (NM_019358) Rat Tagged ORF Clone
Tag: Myc-DDK
Symbol: Pdpm
Synonyms: E11; Gp38; OTS-8; RTI40; T1-alpha
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >RR213710 representing NM_019358
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGTGGACCGCCAGTGTGCTCTGGGTTTTGGGAGCGTTTGGTTCTGGGACTCAGCGCAGGGAGGGG
 CCATCGGTGCGCTAGAAGATGATCTTGTGACCCCGGGTCCAGGAGATGACATGGTGAACCCAGGTCTAGA
 AGACAGAAATAGAAACCACGGACTACAGGAACTGGATAAGTCAACTGCCAAGGCGCCTCTGGTACCA
 ACACAGCCTCCCATTGAGGAACTGCCGACCTCAGGAACCTCCGACCACGATCACAAAGAACATGAGAGTA
 CGACCACTGTCAAAGCGGTGACCAGCCACTCCACGGACAAGAAAACAACCTCACCCAATAGAGATAACGC
 AGGCGGTGAAACGACAGACCAGATAAGAAAGATGGCTTGGCTGTAGTGACCCTGGTCGGAATCATAATT
 GGCGTCTTGTAGCCATTGGCTTCATTGGAGGGATCATCATTGTGGTTATGAGGAAGATTCTGGAAGGT
 TCTCGCCC

ACGCGTACGCGGCGGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RR213710 representing NM_019358
 Red=Cloning site Green=Tags(s)

MWTAPVLLWVLGSVWFWSAQGGAIGALEDLVTGPGDDMVNPGLEDRIETDITGELDKSTAKAPLVP
 TQPIIEELPTSGTSDHDHKEHESTTTVKAVTSHSTDKKTTHPNRDNAGGETQTDDKKDGLAVVTLVGIII
 GVLLAIGFIGGIIIVMRKISGRFSP

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Restriction Sites: SgfI-MluI



Cloning Scheme:



Plasmid Map:



ACCN: NM_019358

ORF Size: 498 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_019358.1</u> , <u>NP_062231.1</u>
RefSeq Size:	1854 bp
RefSeq ORF:	501 bp
Locus ID:	54320
UniProt ID:	<u>Q64294</u>
Cytogenetics:	5q36
MW:	17.6 kDa
Gene Summary:	membrane glycoprotein that may function in lung development [RGD, Feb 2006]