

## Product datasheet for **RR213713**

### Ube2d1 (NM\_001108530) Rat Tagged ORF Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Ube2d1 (NM\_001108530) Rat Tagged ORF Clone  
**Tag:** Myc-DDK  
**Symbol:** Ube2d1  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**ORF Nucleotide Sequence:** >RR213713 representing NM\_001108530  
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**GCGATCGCC**

ATGGCGCTGAAGCGGATCCAGAAAGAATTAAGTGATTTACAGCGTGATCCACCTGCCACTGCTCCGGG  
GACCCGTGGGAGATGACTTGTTCCATTGGCAAGCAACCATCATGGGGCCCCGACAGCGCTATCAAGG  
TGGAGCTTCTTCTCACTGTCCACTTTCCGACAGACTATCCCTTAAACCACCAAGATTGCTTTCACA  
ACAAAAATCTACCACCAAAATATAAACAGTAACGGGAGTATTTGTCTTGACATCCTGAGGTCCCAGTGGT  
CGCCTGCTTTGACCGTATCGAAAGTCTTGTATCAATATGCTCCCTGCTCTGTGACCCTAATCCAGATGA  
TCCCTTAGTACCAGATATTGCACAGATCTATAAATCAGACAAAGAAAAATACAACAGGCACGCGAGAGAA  
TGGACGCAGAAATACGCAATG

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Protein Sequence:** >RR213713 representing NM\_001108530  
Red=Cloning site Green=Tags(s)

MALKRIQKELSDLQRDPPAHCSAGPVGDDLFWQATIMGPPDSAYQGGVFFLTVHFPTDYPFKPPKIAFT  
TKIYHPNINSNGSICLDILRSQWSPALTVSKVLLSICSLLCDPNPDDPLVPDIAQIYKSDKEKYNRHARE  
WTQKYAM

**TR**TRPLEQKLISEEDLAANDILDYKDDDDKV

**Restriction Sites:** SgfI-MluI




[View online »](#)

Cloning Scheme:



Plasmid Map:



ACCN: NM\_001108530

ORF Size: 441 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

<b>OTI Annotation:</b>	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u><a href="#">NM_001108530.1</a></u> , <u><a href="#">NP_001102000.1</a></u>
<b>RefSeq Size:</b>	1360 bp
<b>RefSeq ORF:</b>	444 bp
<b>Locus ID:</b>	361831
<b>UniProt ID:</b>	<u><a href="#">D3ZDK2</a></u>
<b>Cytogenetics:</b>	20p11
<b>MW:</b>	16.6 kDa
<b>Gene Summary:</b>	Accepts ubiquitin from the E1 complex and catalyzes its covalent attachment to other proteins. In vitro catalyzes 'Lys-48'-linked polyubiquitination. Mediates the selective degradation of short-lived and abnormal proteins. Functions in the E6/E6-AP-induced ubiquitination of p53/TP53. Mediates auto-ubiquitination of STUB1, TRAF6 and TRIM63/MURF1. Ubiquitinates STUB1-associated HSP90AB1 in vitro. Lacks inherent specificity for any particular lysine residue of ubiquitin. Essential for viral activation of IRF3. Mediates polyubiquitination of CYP3A4 (By similarity). Mediates ubiquitination of PEX5. [UniProtKB/Swiss-Prot Function]