

Product datasheet for **RR213716L3V**

Camkk2 (NM_031338) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Camkk2 (NM_031338) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Camkk2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_031338
ORF Size:	1761 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR213716).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_031338.1 , NP_112628.1
RefSeq Size:	2132 bp
RefSeq ORF:	1764 bp
Locus ID:	83506
UniProt ID:	O88831
Cytogenetics:	12q16
Gene Summary:	calmodulin-dependent kinase kinase; involved in activating Cam-dependent kinase IV [RGD, Feb 2006]



[View online »](#)