

## Product datasheet for RR213729

### Ccnl2 (NM\_001013094) Rat Tagged ORF Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Ccnl2 (NM\_001013094) Rat Tagged ORF Clone  
**Tag:** Myc-DDK  
**Symbol:** Ccnl2  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**ORF Nucleotide Sequence:** >RR213729 representing NM\_001013094  
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGGCGGGCGGCGGCAGCCGGGGCTTCGGGGTTGATGGCTCCTGCGTTGGCGGCCTGCTCTCCGGCTCGG  
GGGAGCAGCCCCAGGGTCTCAGGGGTGCTGATCGGTGACAGGTATACTCCGGAGTTCTCATCACCTT  
GGAAAAGTGCCTGCTGCCCGACGACAAGCTCCGCTTACGCCATCCATGTCGAGCGGTCTCGACATCGAT  
ACAGAGACTGGACTCCGGTAGTGGGCTGCGAGCTCATCCAGGCAGCCGGCATCCTACTCCGCCTGCCGC  
AGGTGGCCATGGCTACAGGGCAGGTGTTGTTCCAGCGCTTCTTTTACACTAAGTCCTTCTGAAACATTC  
CATGGAGCATGTGTCCATGGCTTGTGTTACCTGGCTTCCAAAATAGAAGAGGCTCCAAGACGAATCAGA  
GATGTCATCAACGTATTCCATCGCCTTCGGCATCTGCGAGAGAAAAAGAAACCTGTGCCTCTGGTCTGG  
ATCAAGAGTATGTTAATCTGAAGAACCAGATTATAAAGGCAGAGAGACGGGTTCTCAAGGAATTGGGCTT  
TTGTGTCCATGTGAAGCATCCTCACAAGATAATCGTTATGTACCTTCAGGTGTTAGAATGTGAGCGTAAT  
CAGCACCTGGTCCAGACTGCATGGGTAGCCTCTGAGGGTAAG

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Protein Sequence:** >RR213729 representing NM\_001013094  
Red=Cloning site Green=Tags(s)

MAAAAAGASGLMAPALAACSSSGSAGAPGSQGVLIIGDRLYSGLVITLENCLLPDDKLRFTPSMSSGLDID  
TETGLRVVGCLEIQAAGILLRLPQYAMATGQVLFQRFFYTKSFVKHSMHVSMACVHLASKIEEAPRRIR  
DVINVFHRLRHLREKKKPVPLVLDQEYVNLKNQIIKAERRVLKELGFCVHVKHPHKIIVMYLQVLECERN  
QHLVQTAWVASEGK


**TR**TRPLEQKLISEEDLAANDILDYKDDDDKV




[View online »](#)

Restriction Sites: SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM\_001013094

ORF Size: 672 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u>NM_001013094.1, NP_001013112.1</u>
<b>RefSeq Size:</b>	1171 bp
<b>RefSeq ORF:</b>	675 bp
<b>Locus ID:</b>	298686
<b>UniProt ID:</b>	<u>Q5I0H5</u>
<b>Cytogenetics:</b>	5q36
<b>MW:</b>	24.6 kDa
<b>Gene Summary:</b>	Involved in pre-mRNA splicing. May induce cell death, possibly by acting on the transcription and RNA processing of apoptosis-related factors.[UniProtKB/Swiss-Prot Function]