

## Product datasheet for **RR213730**

### Srp9 (NM\_001126095) Rat Tagged ORF Clone

#### Product data:

Product Type: Expression Plasmids

Product Name: Srp9 (NM\_001126095) Rat Tagged ORF Clone

Tag: Myc-DDK

Symbol: Srp9

Vector: pCMV6-Entry (PS100001)

E. coli Selection: Kanamycin (25 ug/mL)

Cell Selection: Neomycin

ORF Nucleotide Sequence: >RR213730 representing NM\_001126095  
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGCCTCAGTTCAGACCTGGGAGGAGTTCAGCCGCGCGCTGAGAAGCTCTACTTGGCGGACCCCATGA  
AGGTACGTGTGGTTCTCAAGTACAGGCATGTTGATGGGAATTTGTGTATCAAAGTGACGGACGATCTAGT  
CTGTTTGGTGTACAGAACAGACCAAGCTCAAGATGTAAGAAGATTGAGAAATTCACAGTCAGTTAATG  
CGACTTATGGTGGCCAAGGAATCCCGCAGTGTGCTATGGAAACAGAA

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RR213730 representing NM\_001126095  
Red=Cloning site Green=Tags(s)

MPQFQTWEEFSRAAEKLYLADPMKVRVVLKYRHVDGNLCIKVTDLVCLVYRTDQAQDVKKIEKFHSQLM  
RLMVAKESRSVAMETE

**TR**TRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: Sgfl-MluI




[View online »](#)

Cloning Scheme:



Plasmid Map:



ACCN: NM\_001126095

ORF Size: 258 bp

<b>OTI Disclaimer:</b>	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
<b>OTI Annotation:</b>	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"> <li>1. Centrifuge at 5,000xg for 5min.</li> <li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li> <li>3. Close the tube and incubate for 10 minutes at room temperature.</li> <li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li> <li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li> </ol>
<b>RefSeq:</b>	<a href="#">NM_001126095.1</a> , <a href="#">NP_001119567.1</a>
<b>RefSeq Size:</b>	1328 bp
<b>RefSeq ORF:</b>	261 bp
<b>Locus ID:</b>	690345
<b>Cytogenetics:</b>	13q26
<b>MW:</b>	10.1 kDa