

Product datasheet for **RR213749L3V**

Ube2I6 (NM_001024755) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Ube2I6 (NM_001024755) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Ube2I6
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001024755
ORF Size:	459 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR213749).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001024755.1 , NP_001019926.1
RefSeq Size:	1438 bp
RefSeq ORF:	462 bp
Locus ID:	295704
UniProt ID:	Q4V8J2
Cytogenetics:	3q24
Gene Summary:	Catalyzes the covalent attachment of ubiquitin to other proteins. Functions in the E6/E6-AP-induced ubiquitination of p53/TP53. Promotes ubiquitination and subsequent proteasomal degradation of FLT3.[UniProtKB/Swiss-Prot Function]



[View online »](#)