

Product datasheet for **RR213838**

Pycr1 (NM_001105857) Rat Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Pycr1 (NM_001105857) Rat Tagged ORF Clone
Tag: Myc-DDK
Symbol: Pycr1
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >RR213838 representing NM_001105857
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGAGCGTAGGCTTCATCGGGGCAGGGCAGCTAGCCTTTGCCCTAGCCAAGGGCTTCACAGCTGCAGGTG
TTCTGGCTGCTCACAAGATAATGGCCAGCTCCCCAGACATGGACCAAGCCACAGTCTCTGCCCTCCGGAA
GATAGGGGTGAACCTGACACCCCAACAAGGAGACGGTGTCCACAGCGATGTGCTCTTCCTGGCTGTG
AAGCCCCACATCATCCCCCTTCATCCTGGATGAAATTGGTGTAACTTGAGGACAGGCACATTGTAGTGT
CCTGTGCGGCAGGTGTCACCATCAACTCCATTGAAAAGAACTGACAGCGTTCCAACCTGCTCCCAAAGT
CATCCGATGCATGACTAACACCCCGGTCGTGGTACGGGAGGGGTTACTGTGTACGCCACGGTACTCAC
GCTCAGGTGGAAGATGGCAAGCTTGTGGAGCAGCTCATGGGCAGCGTGGGCTTCTGCACAGAGGTGGAAG
AAGATCTGATCGATGCCGTACGGGGCTCAGTGGCAGCGGCCCTGCCTATGCTTTCACAGCCCTGGATGC
TCTGGCTGATGGTGGTGTAAAAATGGGGCTTCCAAGACGTCTGGCAGTCCGCCTTGGGGCCAGGCCCTC
CTGGGGGCAGCCAAGATGCTACTAGATTGAGAACAACATCCGGGCCAGCTCAAGGACAATGTCTGCTCTC
CTGGTGGGGCCACCATCCATGCCTTGCATGTGCTGGAGAGTGGGGCTTTCGTTCCCTGCTTATCAATGC
TGTGGAAGCCTCCTGTATCCGCACACGGGAGCTGCAGACCATGGCTGATCAGGAAACCGTCTCCCTGCT
GCCATAAAGAAGACTGTACTGGACAAGGTGAAGCTAGACTCCTCTGCTGGGTCTCTGTCTTCTGGCC
ATGTTAAGCCGATGTCC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



Protein Sequence: >RR213838 representing NM_001105857
Red=Cloning site Green=Tags(s)

MSVGFIGAGQLAFALAKGFTAAGVLAHKIMASSPMDQATVSALRKIGVNLTPHNKETVCHSDVFLAV
 KPHIIPFILDEIGANIEDRHIVVSCAAGVTINSIEKKLTAFQPAPKVIKRCMTNTPVVVREGVTVYATGTH
 AQVEDGKLVLEQLMGVSVGFCTEVEEDLIDAVTGLSGSGPAYAFTALDALADGGVKMGLPRRLAVRLGAQAL
 LGAAKMLLDSEQHPGQLKDNVCSPPGATIHAIHLVLESGGFRSLLINAVEASCIRTRELQTMADQETVSPA
 AIKKTVLDDKVKLDSSAGVSLSSGHVKPMS


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_001105857

ORF Size:	927 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001105857.2 , NP_001099327.1
RefSeq Size:	930 bp
RefSeq ORF:	930 bp
Locus ID:	287877
Cytogenetics:	10q32.3
MW:	32.2 kDa