

Product datasheet for **RR213851L4V**

Tceal1 (NM_001009675) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Tceal1 (NM_001009675) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Tceal1
Synonyms:	MGC105656
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001009675
ORF Size:	495 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR213851).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001009675.1 , NP_001009675.1
RefSeq Size:	1215 bp
RefSeq ORF:	498 bp
Locus ID:	302593
UniProt ID:	Q5PPP3
Cytogenetics:	Xq32
Gene Summary:	May be involved in transcriptional regulation. Modulates various viral and cellular promoters in a promoter context-dependent manner. Does not bind DNA directly (By similarity). [UniProtKB/Swiss-Prot Function]



[View online »](#)