

## Product datasheet for **RR213869L3V**

### Trap1 (NM\_001039001) Rat Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Trap1 (NM_001039001) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Trap1
Synonyms:	MGC95278
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001039001
ORF Size:	2118 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR213869).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001039001.1</a> , <a href="#">NP_001034090.1</a>
RefSeq Size:	2264 bp
RefSeq ORF:	2121 bp
Locus ID:	287069
UniProt ID:	<a href="#">Q5XHZ0</a>
Cytogenetics:	10q12



[View online »](#)

**Gene Summary:**

Chaperone that expresses an ATPase activity. Involved in maintaining mitochondrial function and polarization, downstream of PINK1 and mitochondrial complex I. Is a negative regulator of mitochondrial respiration able to modulate the balance between oxidative phosphorylation and aerobic glycolysis. The impact of TRAP1 on mitochondrial respiration is probably mediated by modulation of mitochondrial SRC and inhibition of SDHA.  
[UniProtKB/Swiss-Prot Function]