

Product datasheet for **RR213900**

Spaca5 (NM_001108058) Rat Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Spaca5 (NM_001108058) Rat Tagged ORF Clone
Tag: Myc-DDK
Symbol: Spaca5
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >RR213900 representing NM_001108058
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGCAGTCTGCAGTACTGTGGTGGTATCTTGGCCACACTGTTGCTTCCACAGTGGATGCCAAGATT
ATGAACGCTGTGAGCTGGCAAGGAAGCTGGAGAAAGCTGGCCTCAATGGCTTCAAAGGCTACACTGTTGG
AGACTGGCTGTGTGGCACACTATGAAAGTGGTTTTGACACCTCTTTTGTGGACCACAATCCAGATGGC
AGCAGTGAATATGGCATTTCAGTTGAACTCTGCTTGGTGGTGTAAATGGCATCACACCCACTCAGA
ACCTCTGCCACATGGATTGTAATGACCTACTCAACCGCCATATTCTGGATGATATCATGTGTGCCAAGAG
GGTTGTATCCTCACATAAGAGTATGAAGGCCTGGGATTCTGGACCCGGCACTGTGCTGGTATGATTTA
TCTGAATGGCTCAAGGGGTGTGATGTGCCTCTGAAAGCTGACTCAAACCTATGATGACTGG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RR213900 representing NM_001108058
Red=Cloning site Green=Tags(s)

MQFCSTVVVILATLLLATVDAKIYERCELARKLEKAGLNGFKGYTVGDWLCVAHYESGFDTSFVDHNPDG
SSEYGIFQLNSAWWCNNGITPTQNLCHMDCNDLLNRHILDDIMCAKRVSSSHKSMKAWDSWTRHCAGHDL
SEWLKGC DVPLKADSNYDDW

TRTRPLEQKLI SEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI




[View online »](#)

Cloning Scheme:



Plasmid Map:



ACCN: NM_001108058

ORF Size: 480 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001108058.1, NP_001101528.1</u>
RefSeq Size:	553 bp
RefSeq ORF:	483 bp
Locus ID:	314431
Cytogenetics:	Xq11
MW:	18.1 kDa