

Product datasheet for **RR213902**

LOC361016 (NM_001014137) Rat Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	LOC361016 (NM_001014137) Rat Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	LOC361016
Synonyms:	MGC116196
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

**ORF Nucleotide
Sequence:**

>RR213902 representing NM_001014137
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGC**

ATGTTCAAAAATGAGTACCAGGGAGGCGCATTTGTTGAAATTTTTCAGCGCTCAGGGGAAAAATCCTGGAG
 CAAAGTGAAGATCTTGGGCAGTCCATCTGTGATCTGGAAAGAGTTTGATAAAGAAGTTAAAGCTTTGT
 GTTTGCTCTGGAAGGCAGCAGCCAAACACAGAATCCAGTTGCCGAAGGAGAATAAGCAAAATCCTTGGA
 CTGATCCAGAGGTTTCTGTGCTACAGATTTACGTGCCCTGGGACAGGACTTCTCCACCGAATTACTAA
 TTAAGTACTTAGGAAACATCAAAAGAAGTTGTATTTATCAACAGTCCATAAGGAAGTGTCTTCAACACC
 TCTTCACGCAAAAATCCCACTCTTCATGATCAAGCGTAAATCTGGTGCAATCTGTGCATTGACCTGGTC
 GCCTTCACCAAGTGAATATTCAAAGGAGCAGTTTTCCAGTCGCTGGATGGAATTATCGTCTCTGCTAACT
 GTAAGCTGCGGAAGATCTTCACTTTGAAATTCAAGCCTCGAGAGACGGCCGATAGAGACGATGAACCCGC
 AGATATTATTCCACGAAGCTGTCAACTCGCAACAGATGTTCCGCTGTCACACAGCTGCTTAACATGACA
 AAACCTTGCCAAACTGAAATTAATTTGGAGGTCACTCTAAGATCTGCAGAATCAGACCAGTTTCATCA
 GCAGCAGAGCAGGCAAGTGTCCGGAGCAGTAAAAGTCAAGATGTCTGCCACATTGCCTTTGGGTCCAGGGT
 TCTTGGACCACCTCCACCTCTGGCAGACGGAATAACCTAAGGCTTTCCGCTGAGACAGTGAATCCATT
 GGGTACAAAAATAACCAAGTCTGTCAGCAGCCTGCGGAAGAGAAGCATGTTAACAGTGCAGACACATCAG
 CCTTGCTGATACCTCTGTCTGAACAGCAAGGAGAGAAAGGAAGTTCTCACCCCGTGAAGCAGACACGCC
 TCTACCTGCCAGTCTCCCTACTCCATATCCACATCCCCCTCGAGACCCGTGAGCAGATAAGGGTAGCAGC
 AGAAGAGGACTTGGCCTAAGGAGCAGCAGTGAAGCAGGACAGAGGCGCGCTGTGGGAGTTCTCAGGAA
 ATAGCAGGAGTGAAGAAGACGTAGCCATGGCCGCAACAATCTCACACAGCTGACGACAGAGCCAGTGAG
 AGCCAGCACTCCAGAACCCCGGCTGCTGAACCAACAGATGAGTGGATATTCCCTGAAAGCAGTGGCCCA
 CCCTCCGGCGCAGCAGGTAGCAGCTCATCTTTGCTCCTAGATGGGACTCCCGTACCGCATCACACCTGT
 GGCTCGAAACCAGCAAGGACAGTGACCTGACCAAGGAGAGGGAACCCAGATGGCTCCCAAGGACATTTT
 CACCTTTTCGTCCAGACCAAGATCAGCGCCTCATGGGAAGAGTACGATCTGTCCCCACGGGTGCCCA
 TGTATTCTGGACCCAAAAGAGGACAGCCGAGTGACCAAGGAGACACCGAGTTGGAAGACGATTTTATG
 GCAGTGACAGCAGCAAGAGGAATACAAGTGGAGAACTATCAGCCAGCCAGATGAGCGAATCTGAGCT
 ACAGATGCTGGCAAGCCTCAGGCGGCAGCAGAATGAAGACCTAGAGGACACCGGGGCTCCTCATGGCCTG
 AGCGCATCCAGGTGGACAAGTGAATGTCAGCATCAGTACCAGCAGCGATGACACGACCAGTGGAAGT
 CCTGTCTCCACCTCCTGTGAACCAAGGCGCTCACTATCAGAAAGAAATGAACCCCTTCTCCATCTAA
 CCCTCGAGACTGGTTGAACATGTTGAGCCACCAATTGTCCCTCCAGTACGAGCCGCTGGAGCAGTCT
 CTGGACTCCTCTGCAAGTTTGAGTGTGCAAGGAGGAGCTCTGCTACTGACTCAGCCCCACCTCGTTCTCT
 GTGCAGTTGGGTGT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RR213902 representing NM_001014137
 Red=Cloning site Green=Tags(s)


MFKNEYQGGAFFEISAQGKNPGAKWKILGSPSVIWKEFDKEVKSFFVFLGSSQTNRIQLPKENKQILG
 LIQRFLVLQIYVPLGQDFSTELLITDLGNIKRRLYLSTVHKEVSSTPLHAKIPLFMKRKIWCNLCIDL
 AFTSEIFKGAVFQSLDGIIVSANCKLRKIFTLKFKPRETADRDDEPADIIIPRSCQLATDVPHVTQLLNM
 KLRQTEIKFGGHLPLRSAESDQFISSRAGSVRSSKSQDVCHIAFGSRVLGPPPPSGRRNNLRLSAETVRSI
 GYKNNQSCQQPAEEKHVNSADTSALLIPLSEQQGEKGSSHPVKQTTPLPASLPTYPHPPRDPADKGS
 RRGLGLRSSSGSRTEARCGSSSGNSRSEEDVAMAANNLTQLTTEPVRASTPEPPAAEPPDEWIFPESSGP
 PSGAAGSSSSLLLDGDSRTASHLWLETSKSDSDPDQAEQTMAPKDIPTFSSRPRSAPHGKSHDLSPGCP
 CILDPKEDSRVTKGDTLEDDFYGSDSSEEEYNWRNYQPSQMSESELQMLASLRRQNEDLEDTGAPHGL
 SASQVDNCNVSISTSSDDTTTWNCLPPPVNQGRHYQKEMNPPSPSNPRDWNMLSPPIVPPSQPLEQS
 LDSSASLSVQGGALLLTQPHLSCAVGC

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites:

SgfI-MluI

Cloning Scheme:



ACCN: NM_001014137

ORF Size: 1974 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001014137.2](#), [NP_001014159.2](#)

RefSeq Size: 2306 bp

RefSeq ORF: 1977 bp


Locus ID: 361016

UniProt ID: [Q4V7B1](#)

Cytogenetics: 15p14

MW: 71.9 kDa

Product images:



Circular map for RR213902