

Product datasheet for **RR213925**

Rxfp3 (NM_001008310) Rat Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Tag:	Myc-DDK
Symbol:	Rxfp3
Synonyms:	GPCR135; RGD1311550; Rln3r1; Salpr
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide Sequence: >RR213925 representing NM_001008310
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCCGCGATCGCC

ATGCCAAAGCGCACCTGAGCATGCAAGTGGTCTTCAACACCGCAGCCCCATGAGTAAGGCAGCTG
 CGGGTGATGAGCTCTCCGATTCTTCGGCTGATCCCAGACTTGCTGGAGGTTGCCAACAGGAGCAGCAA
 TGCCTGCTGACGCTTCAAGACTTGTGGTGGGAGCTGGGGCTGGAGTTGCCCGACGGTGCGCGCCTGGG
 CATCCCCGGGCGAGCGTGGGGCAGAGAGCGCGGACACAGAGGCCAGGGTACGGATCCTCATCAGCGCCG
 TTTACTGGGTGGTTGTGCCCTGGGACTGGCTGGCAACCTGCTGGTTCTCTACCTGATGAAGAGCAAACA
 GGGCTGGCGCAAACTCTCCATTAACTCTTTGTCACTAACCTGGCGCTGACTGACTTTTCAGTTTGTGCTC
 ACTCTGCCCTTCTGGGCGGTGGAGAACGCACTAGATTTCAGTGGCCCTTTGGCAAGGCCATGTGTAAGA
 TCGTATCTATGGTGACATCCATGAACATGTATGCCAGCGTCTTCTTTCTACTGCTATGAGTGTGGCGCG
 CTACCACTCGGTGGCCTCAGCTCTCAAGAGCCATCGGACCCGCGGCATGGCCGTGGCGACTGCTGCGGC
 CAGAGCTTGGGGGAGAGCTGCTGTTTCTCAGCCAAGGTGCTGTGTGGATTGATCTGGGCTTCTGCCGCGA
 TAGCTTCGCTGCCCAATGTCAATTTTTCTACCACCATCAATGTGTTGGGCGAGGAGCTGTGCCTCATGCA
 CTTTCCGGACAAGCTCCTGGGTTGGGACCGGCAGTTCTGGCTGGGTTGTACCACCTGCAGAAGGTGCTG
 CTGGGCTTCTGCTGCCGCTGAGCATCATAGTTTGTGTTACCTGTTGCTCGTGCCTTCTCTCCGACC
 GCCGCGTAGTGGGGACAACGGATGGAGCAACAGCGCCTGGGGGGAGCCTGAGTACAGCCGGCGCTCGGAG
 ACGCTCCAAGGTCACCAAGTGGTGACCATCGTAGTCCTTTCTTCTTATGTTGGCTGCCCAACCAA
 GCGCTACCACTGGAGCATCCTCATAGTTCAACGTAGTCCCTTCAGTCAGGAGTACTTTTCAGTGCC
 AAGTGACGCTTCCAGTCAGCGTGTGCTGGCACACTCCAACAGCTGCCTCAACCCCATCCTCTACTG
 CTTAGTGCGCCGCGAGTTCCGCAAGGCGCTCAAGAACCTGCTGTGGCGTATAGCATCGCCTTCGCTCACC
 AGCATGCGCCCTTACCGCCACCACCAAGCCAGAACCTGAAGATCACGGGCTGCAGGCCCTGGCGCCAC
 TTAATGCTACTGCAGAGCTGACCTGATCTACTATCCACCCGGTGGTGGTCTACAGCGGAGGTGCTA
 CGACCTTCTCCCTAGCAGCTCTGCCTACCGACCTTCTCCCTAGCAGCTCTGCCTAC

AGCGGACCGACGCGTACGCGCGGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCC
 TGGATTACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RR213925 representing NM_001008310
 Red=Cloning site Green=Tags(s)


MPKAHLSMQVASATTAAPMSKAAAGDELSGFFGLIPDLLLEVANRSSNASLQLQDLWWELGLELPDGAAPG
 HPPGSGGAESADTEARVRILISAVYVVCALGLAGNLLVLYLMSKQGWKSSINLFVTNLALTDQFVL
 TLPFWAVENALDFKWPFGKAMCKIVSMVTSMNMYASVFFLTAMSVARYHSVASALKSHRTRGHGRDCCG
 QSLGESCCFSKVLCLIWASAAIASLPNIVFSTTINVLEELCLMHFPDKLLGWDRQFWLGLYHLQKVL
 LGFLLPLSIISLCYLLVRFISDRRVVGTDTGATAPGGSSTAGARRRSKVTKSVTIVVLSFFLCWLPNQ
 ALTTWSILIKFNVPFVSQEFYQCQVYAFVSVCLAHSNSCLNPILYCLVRREFRKALKNLLWRIASPSLT
 SMRPFTATTKEPEDHGLQALAPLNATAEPDLIYPPGVVVYSGGRYDLLPSSSAYRPS*QLCL

SGPTRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/ja1686_g07.zip

Restriction Sites: SgfI-RsrII

Cloning Scheme:



ACCN: NM_001008310

ORF Size: 1456 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq Size: 1431 bp


RefSeq ORF: 1431 bp

Locus ID: 294807

Cytogenetics: 2q16

MW: 52 kDa

Product images:



Circular map for RR213925