

Product datasheet for **RR213930**

Pusl1 (NM_001108699) Rat Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Pusl1 (NM_001108699) Rat Tagged ORF Clone
Tag: Myc-DDK
Symbol: Pusl1
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
ORF Nucleotide Sequence: >RR213930 representing NM_001108699
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGCAGGAGGCTGCCAGCATCTCCTTGGGACACATGATTTTGTGCCTTCCAGTCTGCTGGCAGCCAG
 TCCACATGCTGTACGCACACTGCGAAGAGTCTCTGTGTCTCCTGGTCCAGCCAGCTTATTTGTCTCCC
 TCAGGAGAGCAGGAGGCTGCAGTTCTGGACCTTGGAGTTTGAAGCCAGTCTTTCCTCTATCGACAGGTG
 CGGAGGATGACAGCTGTGCTGGTGGCTGTGGGGCTAGGGATTCTGGCACCTACACAAGTAAAGTGATTC
 TGGAGAGTCAAGATCCCTTGGGCAAGTATCAGGCTCGAGTTGCCCCAGCCATGGCCTGTTCTTGAAGTC
 AGTGCTGTATGACAACTTGGCTGGGAAAAC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RR213930 representing NM_001108699
 Red=Cloning site Green=Tags(s)


MQEAAQHLLGTHDFSAFQSAGSPVPHAVRTLRRVSVSPGPASLFVLPQESRRLQFWTLEFESQSFLYRQV
 RRMTAVLVAVGLGILAPTQVKVILESQDPLGKYQARVAPAHGLFLKSVLYDNFGWEN

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI


[View online »](#)

Cloning Scheme:



ACCN: NM_001108699

ORF Size: 381 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method: 1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.


RefSeq: [NM_001108699.1](#), [NP_001102169.1](#)

RefSeq Size: 1688 bp

RefSeq ORF: 384 bp

Locus ID: 362681
 Cytogenetics: 5q36
 MW: 14.2 kDa

Product images:



Circular map for RR213930